








- All patient handling equipment should be inspected prior to use. Follow manufacturer's guidelines for care, cleaning and inspection.

EQUIPMENT	TYPE OF TRANSFER	PATIENT MOBILITY ASSESSMENT	SLING TYPE
<p><u>Transfer Belt</u></p> 	<ul style="list-style-type: none"> • Sit to Sit 	Minimal Assist	N/A
<p>Sara Steady</p> 	<ul style="list-style-type: none"> • Sit to Sit 	Patient able to pull themselves up, not able to step/walk.	N/A
<p><u>Sit/Stand Lift</u></p> 	<ul style="list-style-type: none"> • Sit to Sit 	Use if patient partial strength in at least one leg and one arm, (two legs if bariatric)	Harness
<p><u>Floor Based Total Lift</u></p> 	<ul style="list-style-type: none"> • Bed/chair/commode etc. • Lateral Transfer • Wound Care • Patient Turns 	Maximal Assist or Uncooperative	<ul style="list-style-type: none"> • Transfer Sling • Walking Harness • Hygiene Sling • Repositioning Sling • Limb Sling • <u>Sling Types</u> • <u>General Guidelines For Sling Fitting</u>
<p>Bariatric Floor Based Total Lift</p> 	<ul style="list-style-type: none"> • Bed/chair/commode etc. • Lateral Transfer • Wound Care • Patient Turns 	Maximal Assist or Uncooperative	<ul style="list-style-type: none"> • Sized Transfer Sling • Sized Repositioning Sling • Limb Sling • <u>Sling Types</u> • <u>General Guidelines For Sling Fitting</u>

<p><u>Ceiling Lift</u></p> 	<ul style="list-style-type: none"> • Bed/chair/commode etc. • Lateral Transfer • Wound Care • Patient Turns 	<p>Maximal Assist or Uncooperative</p>	<ul style="list-style-type: none"> • Transfer Sling • Hygiene Sling • Repositioning Sling • Limb Sling <ul style="list-style-type: none"> • Sling Types • General Guidelines For Sling Fitting
<p>Bariatric Twin Module Ceiling Lift</p>	<ul style="list-style-type: none"> • Bed/chair/commode etc. • Lateral Transfer • Wound Care • Patient Turns 	<p>Maximal Assist or Uncooperative</p>	<ul style="list-style-type: none"> • Transfer Sling • Repositioning Sling • Limb Sling
<p>Air-Assisted Device (HoverJack)</p> 	<ul style="list-style-type: none"> • Recovery from floor • Lateral Transfer • Hovermatt/jack SWP 	<p>Maximal Assist</p>	<p>N/A</p>

Links:

- [WorkSafe BC Bulletin – Ceiling Lifts vs Floor Lifts](#)
- [WorkSafe BC Bulletin – Ceiling Lifts](#)